




UIF rabbit pAb

Catalog No	YP-Ab-11940
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	FYTDD1 UIF
Immunogen	Synthesized peptide derived from human UIF AA range: 50-100
Specificity	This antibody detects endogenous levels of UIF at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Forty-two-three domain-containing protein 1; Protein 40-2-3; UAP56-interacting factor; UIF; FYTDD1 UIF
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus, nucleoplasm . Nucleus speckle .
Tissue Specificity	Expressed in a wide variety of cancer types.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night