




FA69C rabbit pAb

Catalog No	YP-Ab-11941
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM69C C18orf51
Immunogen	Synthesized peptide derived from human FA69C AA range: 345-395
Specificity	This antibody detects endogenous levels of FA69C at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C18orf51;DIPK1C;Divergent protein kinase domain 1C;Protein FAM69C;FA69C; FAM69C C18orf51
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum membrane ; Single-pass type II membrane protein .
Tissue Specificity	
Function	
Background	This gene encodes a member of the FAM69 family of cysteine-rich type II transmembrane proteins. These proteins localize to the endoplasmic reticulum but their specific functions are unknown. [provided by RefSeq, Nov 2011],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night