




COLL2 rabbit pAb

Catalog No	YP-Ab-11942
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	IHC; IF
Gene Name	CLPSL2 C6orf126 UNQ3045/PRO9861
Immunogen	Synthesized peptide derived from human COLL2 AA range: 37-87
Specificity	This antibody detects endogenous levels of COLL2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	COLL2; CLPSL2 C6orf126 UNQ3045/PRO9861
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity: Belongs to the colipase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).