




RTKN2 rabbit pAb

Catalog No	YP-Ab-11944
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	RTKN2 PLEKHK1
Immunogen	Synthesized peptide derived from human RTKN2 AA range: 473-523
Specificity	This antibody detects endogenous levels of RTKN2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PH domain-containing family K member 1;Pleckstrin homology domain-containing family K;RTKN2;RTKN2 PLEKHK1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	Expressed in lymphocytes, CD4 positive T-cells and bone marrow-derived cells. Also expressed in lung, colon, thymus and brain.
Function	function:May play an important role in lymphopoiesis.,induction:Up-regulated in cells resistant to 25-hydroxy cholesterol (25-OHC). Down-regulated in lymphocytes activated by treatment with phytohemagglutinin (PHA)., similarity:Contains 1 PH domain.,tissue specificity:Expressed in lymphocytes, CD4 positive T-cells and bone marrow-derived cells. Also expressed in lung, colon, thymus and brain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HEK293 cells, primary antibody was diluted at 1:1000, 4° over night