




F102A rabbit pAb

Catalog No	YP-Ab-11949
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM102A C9orf132 EEIG1
Immunogen	Synthesized peptide derived from human F102A AA range: 72-122
Specificity	This antibody detects endogenous levels of F102A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FAM12A;F102A;FAM102A C9orf132 EEIG1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	function:May play a role in estrogen action.,induction:By 17-beta-estradiol but also by a group of natural and synthetic estrogens as well as by estrogenic environmental compounds. Repressed by the antiestrogen 4-hydroxy-tamoxifen., similarity:Belongs to the FAM102 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4° over night