



SAHH3 rabbit pAb


Catalog No	YP-Ab-11969
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	AHCYL2 KIAA0828
Immunogen	Synthesized peptide derived from human SAHH3 AA range: 157-207
Specificity	This antibody detects endogenous levels of SAHH3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AdoHcyase 3;AHCYL2;FLJ21719;KIAA0828;SAHH3;AHCYL2 KIAA0828
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Microsome . Associates with membranes when phosphorylated, probably through interaction with ITPR1. .
Tissue Specificity	
Function	catalytic activity:S-adenosyl-L-homocysteine + H(2)O = L-homocysteine + adenosine.,cofactor:Binds 1 NAD per subunit.,cofactor:NAD.,pathway:Amino-acid biosynthesis; homocysteine biosynthesis; L-homocysteine from S-adenosyl-L-homocysteine: step 1/1.,similarity:Belongs to the adenosylhomocysteinase family.,
Background	The protein encoded by this gene acts as a homotetramer and may be involved in the conversion of S-adenosyl-L-homocysteine to L-homocysteine and adenosine. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2011],

matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night