




KDEL2 rabbit pAb


Catalog No	YP-Ab-11971
Isotype	IgG
Reactivity	Human;Rat
Applications	WB;ELISA;IHC
Gene Name	KDEL2 UNQ1904/PRO4350
Immunogen	Synthesized peptide derived from human KDEL2 AA range: 414-464
Specificity	This antibody detects endogenous levels of KDEL2 at Human/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Protein O-xylosyltransferase POGLUT3;Protein O-glucosyltransferase 3; POGLUT3;KDEL motif-containing protein 2;EC:2.4.2.-;KDEL2;KDEL2 UNQ1904/PRO4350
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Endoplasmic reticulum lumen .
Tissue Specificity	
Function	similarity:Belongs to the KDEL2 family.,similarity:Contains 1 filamin repeat.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).