




CX059 rabbit pAb

Catalog No	YP-Ab-11972
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	CXorf59
Immunogen	Synthesized peptide derived from human CX059 AA range: 306-356
Specificity	This antibody detects endogenous levels of CX059 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CXorf59;Cilia and flagella associated protein 47;CX059
Observed Band	
Calculated Molecular Weight	
Cell Pathway	integral component of membrane,
Tissue Specificity	
Function	
Background	While this gene is well-supported by transcript data, no functional information on its protein product is currently available. [provided by RefSeq, Dec 2009],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4° over night