




NOP9 rabbit pAb

Catalog No	YP-Ab-11983
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NOP9 C14orf21 KIAA2021
Immunogen	Synthesized peptide derived from human NOP9 AA range: 486-536
Specificity	This antibody detects endogenous levels of NOP9 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	chromosome 14 open reading frame 21;hypothetical protein LOC161424; KIAA2021;NOP9;NOP9 nucleolar protein;NOP9 nucleolar protein homolog; nucleolar protein 9;Pumilio domain containing protein C14orf21;NOP9 C14orf21 KIAA2021
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 6 pumilio repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night