




DSEL rabbit pAb

Catalog No	YP-Ab-12002
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DSEL C18orf4 NCAG1
Immunogen	Synthesized peptide derived from human DSEL AA range: 676-726
Specificity	This antibody detects endogenous levels of DSEL at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	DSEL;DSEL C18orf4 NCAG1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Multi-pass membrane protein .
Tissue Specificity	Expressed in different brain areas as well as in multiple other peripheral tissues.
Function	similarity:Belongs to the dermatan-sulfate isomerase family.,tissue specificity:Expressed in different brain areas as well as in multiple other peripheral tissues.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4° over night