




TVC rabbit pAb

Catalog No	YP-Ab-12005
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	TRGV3 TCRB
Immunogen	Synthesized peptide derived from human TVC AA range: 79-129
Specificity	This antibody detects endogenous levels of TVC at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TVC;TRGV3 TCRB
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane .
Tissue Specificity	
Function	miscellaneous:This is an aberrantly rearranged nucleotide sequence. Two bases are deleted at nucleotide position 357 in this clone. However, the protein sequence can be deduced from the mRNA sequence.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from THP-1 cells,
primary antibody was diluted at 1:1000, 4° over night