




MET14 rabbit pAb


Catalog No	YP-Ab-12013
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;ELISA;IHC
Gene Name	METTL14 KIAA1627
Immunogen	Synthesized peptide derived from human MET14 AA range: 360-410
Specificity	This antibody detects endogenous levels of MET14 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	hMETTL14;Methyltransferase-like protein 14;METT L14;N6-adenosine-methyltransferase non-catalytic subunit;MET14;METTL14 KIAA1627
Observed Band	
Calculated Molecular Weight	50kD
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:Probable methyltransferase.,similarity:Belongs to the MT-A70-like family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).