



TRI47 rabbit pAb

Catalog No	YP-Ab-12016
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TRIM47 GOA RNF100
Immunogen	Synthesized peptide derived from human TRI47 AA range: 199-249
Specificity	This antibody detects endogenous levels of TRI47 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	E3 ubiquitin-protein ligase TRIM47;EC:2.3.2.27;GOA;RNF100;tripartite motif containing 47;TRI47;TRIM47 GOA RNF100
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Nucleus .
Tissue Specificity	Low expression in most tissues. Higher expression in kidney tubular cells. Overexpressed in astrocytoma tumor cells.
Function	similarity:Belongs to the TRIM/RBCC family.,similarity:Contains 1 B box-type zinc finger.,similarity:Contains 1 B30.2/SPRY domain.,similarity:Contains 1 RING-type zinc finger.,tissue specificity:Low expression in most tissues. Higher expression in kidney tubular cells. Overexpressed in astrocytoma tumor cells.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images