




F132A rabbit pAb

Catalog No	YP-Ab-12020
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM132A C1QDC2
Immunogen	Synthesized peptide derived from human F132A AA range: 193-243
Specificity	This antibody detects endogenous levels of F132A at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	F132A;FAM132A C1QDC2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	[Adipolin fC1QTNF12]: Secreted .; [Adipolin gC1QTNF12]: Secreted . In sera is the predominant form. .
Tissue Specificity	Predominantly expressed by adipose tissues.
Function	similarity:Belongs to the FAM132 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night