




TB15A rabbit pAb

Catalog No	YP-Ab-12023
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	IHC;IF
Gene Name	TMSB15A TMSL8 TMSNB
Immunogen	Synthesized peptide derived from human TB15A AA range: 1-24
Specificity	This antibody detects endogenous levels of TB15A at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TB15A;TMSB15A TMSL8 TMSNB
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton .
Tissue Specificity	Neuroblastoma-specific.
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).