




FKB15 rabbit pAb

Catalog No	YP-Ab-12024
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FKBP15 KIAA0674
Immunogen	Synthesized peptide derived from human FKB15 AA range: 847-897
Specificity	This antibody detects endogenous levels of FKB15 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	133 kDa FK506 binding protein;133 kDa FKBP;FK506 binding protein 15;FKBP 133;FKBP 15;FKB15;FKBP15 KIAA0674
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm . Cell projection, axon . Early endosome . Present in axons and neuronal growth cones. .
Tissue Specificity	
Function	domain:The PPlase FKBP-type domain seems to be inactive both for FK506-binding and enzymatic activity.,function:May be involved in the cytoskeletal organization of neuronal growth cones. Seems to be inactive as a PPlase.,similarity:Belongs to the FKBP-type PPlase family.,similarity:Contains 1 PPlase FKBP-type domain.,subcellular location:Present in axons and neuronal growth cones.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night