




ENHO rabbit pAb

Catalog No	YP-Ab-12030
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	ENHO C9orf165 UNQ470/PRO830
Immunogen	Synthesized peptide derived from human ENHO AA range: 15-65
Specificity	This antibody detects endogenous levels of ENHO at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Adropin;C9orf165;energy homeostasis associated;ENHO;UNQ47;ENHO C9orf165 UNQ470/PRO830
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	Expressed in liver and brain.
Function	function:Involved in the regulation of glucose homestasis and lipid metabolism., tissue specificity:Expressed in liver and brain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 30min).