




SERF2 rabbit pAb

Catalog No	YP-Ab-12031
Isotype	IgG
Reactivity	Human; Mouse
Applications	IHC;IF
Gene Name	SERF2 FAM2C
Immunogen	Synthesized peptide derived from human SERF2 AA range: 1-80
Specificity	This antibody detects endogenous levels of SERF2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	IHC-p 1:50-200. IF 1:50-200
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	small EDRK rich factor 2;Protein 4F5-related;Protein 4F5 related;h4F5rel;Gastric cancer-related;SERF2;SERF2 FAM2C
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,cytosol,
Tissue Specificity	
Function	sequence caution:Wrong choice of frame.,similarity:Belongs to the SERF family., subunit:May interact with the N-terminus of HD.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200(room temperature, 45min).