




TM59L rabbit pAb

Catalog No	YP-Ab-12032
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	TMEM59L BSMAP C19orf4
Immunogen	Synthesized peptide derived from human TM59L AA range: 177-227
Specificity	This antibody detects endogenous levels of TM59L at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM59L; TMEM59L BSMAP C19orf4
Calculated Molecular Weight	38kD
Cell Pathway	Golgi apparatus membrane ; Single-pass type I membrane protein .
Tissue Specificity	Expressed preferentially at high level in the brain.
Function	function: May have a role in brain function., similarity: Belongs to the TMEM59 family., subcellular location: Could be associated with an intracellular organelle or membrane., tissue specificity: Expressed preferentially at high level in the brain.,
Background	This gene encodes a predicted type-I membrane glycoprotein. The encoded protein may play a role in functioning of the central nervous system. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night