




KLDC9 rabbit pAb

Catalog No	YP-Ab-12065
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	KLHDC9 KARCA1
Immunogen	Synthesized peptide derived from human KLDC9 AA range: 32-82
Specificity	This antibody detects endogenous levels of KLDC9 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KLDC9;KLHDC9 KARCA1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Contains 3 Kelch repeats.,subunit:Interacts with CCNA1.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night