




WDR76 rabbit pAb

Catalog No	YP-Ab-12071
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	WDR76
Immunogen	Synthesized peptide derived from human WDR76 AA range: 182-232
Specificity	This antibody detects endogenous levels of WDR76 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CDW14;FLJ12973;WD repeat domain 76;WD repeat-containing protein 76;Wdr76;WDR76_HUMAN
Observed Band	
Calculated Molecular Weight	
Cell Pathway	nucleus,
Tissue Specificity	
Function	similarity:Belongs to the WD repeat DDB2/WDR76 family.,similarity:Contains 7 WD repeats.,subunit:Interacts with CUL4A and/or CUL4B.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4°over night