




WDR47 rabbit pAb

Catalog No	YP-Ab-12093
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	WDR47 KIAA0893
Immunogen	Synthesized peptide derived from human WDR47 AA range: 474-524
Specificity	This antibody detects endogenous levels of WDR47 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ90135;KIAA0893;WD repeat domain 47;WDR47;WDR47 KIAA0893
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton . Localization along microtubules is mediated by MAP1S. .
Tissue Specificity	
Function	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data., similarity:Contains 1 CTLH domain.,similarity:Contains 1 LisH domain., similarity:Contains 7 WD repeats.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night