



ATAD5 rabbit pAb

Catalog No	YP-Ab-12109
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ATAD5 C17orf41 FRAG1
Immunogen	Synthesized peptide derived from human ATAD5 AA range: 1054-1104
Specificity	This antibody detects endogenous levels of ATAD5 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ATAD5;ATAD5 C17orf41 FRAG1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:Involved in DNA damage response. Involved in a RAD9A-related damage checkpoint, a pathway that is important in determining whether DNA damage is compatible with cell survival or whether it requires cell elimination by apoptosis. Modulates the RAD9A interaction with BCL2 and thereby induces DNA damages-induced apoptosis.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR. ATR may stimulate the RAD9A dissociation., sequence caution:Contaminating sequence. Potential poly-A sequence., similarity:Belongs to the AAA ATPase family.,subunit:Interacts with RB1 predominantly in G1 phase via its LXCXE motif. Interacts with RAD9A in growing cells. The interaction with RAD9A is reduced after exposure to DNA replication-inhibiting agents.,

Background


matters needing attention

Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HeLa cells, primary antibody was diluted at 1:1000, 4° over night