




PP2D1 rabbit pAb

Catalog No	YP-Ab-12131
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PP2D1 C3orf48
Immunogen	Synthesized peptide derived from human PP2D1 AA range: 9-59
Specificity	This antibody detects endogenous levels of PP2D1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	PP2D1;PP2D1 C3orf48
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night