




WDR53 rabbit pAb

Catalog No	YP-Ab-12139
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	WDR53
Immunogen	Synthesized peptide derived from human WDR53 AA range: 19-69
Specificity	This antibody detects endogenous levels of WDR53 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	WDR53
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the WD repeat WDR53 family.,similarity:Contains 5 WD repeats.,
Background	This gene encodes a protein containing WD domains. The function of this gene is unknown. [provided by RefSeq, Jul 2008],
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night