




CPPED rabbit pAb

Catalog No	YP-Ab-12144
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	CPPED1 CSTP1
Immunogen	Synthesized peptide derived from human CPPED AA range: 68-118
Specificity	This antibody detects endogenous levels of CPPED at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CPPED1; CSTP1; CPPED; CPPED1 CSTP1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm .
Tissue Specificity	Expressed in subcutaneous adipose tissue.
Function	cofactor: Binds 2 divalent metal cations., induction: Transactivated by the large envelope protein of the hepatitis B virus (HBV)., similarity: Belongs to the metallophosphoesterase superfamily. CPPED1 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from SW480 cells,
primary antibody was diluted at 1:1000, 4° over night