




ZBT46 rabbit pAb

Catalog No	YP-Ab-12145
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZBTB46 BTBD4 ZNF340
Immunogen	Synthesized peptide derived from human ZBT46 AA range: 108-158
Specificity	This antibody detects endogenous levels of ZBT46 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Zinc finger and BTB domain-containing protein 46;RINZF;BZEL;BTB-ZF protein expressed in effector lymphocytes;BTBD4;ZBT46;ZBTB46 BTBD4 ZNF340
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May be involved in transcriptional regulation.,similarity:Contains 1 BTB (POZ) domain.,similarity:Contains 2 C2H2-type zinc fingers.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night