




GGNB2 rabbit pAb

Catalog No	YP-Ab-12151
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	GGNBP2 LCRG1 LZK1 ZNF403
Immunogen	Synthesized peptide derived from human GGNB2 AA range: 225-275
Specificity	This antibody detects endogenous levels of GGNB2 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GGNB2; GGNBP2 LCRG1 LZK1 ZNF403
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasmic vesicle . Associated with vesicular structures. .
Tissue Specificity	Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas. Expressed more abundantly in heart, pancreas and skeletal muscle.
Function	function: May be involved in spermatogenesis., miscellaneous: Strongly down-regulated in 40% of primary laryngeal carcinoma and in 6 of 10 various cancer cell lines., subcellular location: Associated with vesicular structures., subunit: Interacts with GGN., tissue specificity: Expressed in heart, brain, placenta, lung, liver, skeletal muscle, kidney and pancreas. Expressed more abundantly in heart, pancreas and skeletal muscle.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

**Usage suggestions**

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from HpeG2 cells, primary antibody was diluted at 1:1000, 4° over night