




S13A1 rabbit pAb

Catalog No	YP-Ab-12158
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	SLC13A1 NAS1 NAS1
Immunogen	Synthesized peptide derived from human S13A1 AA range: 490-540
Specificity	This antibody detects endogenous levels of S13A1 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	SLC13A1; NAS1; hNaSi-1; S13A1; SLC13A1 NAS1 NAS1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane; Multi-pass membrane protein.
Tissue Specificity	Highly expressed in kidney; not detectable in the other tissues tested.
Function	function: Sodium/sulfate cotransporter that mediates sulfate reabsorption in the kidney., induction: Inhibited by thiosulfate, selenate, molybdate, tungstate, citrate and succinate., similarity: Belongs to the SLC13A transporter (TC 2.A.47) family. NADC subfamily., tissue specificity: Highly expressed in kidney; not detectable in the other tissues tested.,
Background	The protein encoded by this gene is an apical membrane Na(+)-sulfate cotransporter involved in sulfate homeostasis in the kidney. Defects in this gene lead to many pathophysiologic problems. [provided by RefSeq, May 2016],
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from A431 cells, primary antibody was diluted at 1:1000, 4° over night