




MILK2 rabbit pAb


Catalog No	YP-Ab-12161
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB;ELISA;IHC
Gene Name	MICALL2
Immunogen	Synthesized peptide derived from human MILK2 AA range: 269-319
Specificity	This antibody detects endogenous levels of MILK2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	JRAB;MICAL L2;MICAL like 2;MICAL like protein 2;MICALL2;MILK2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cell membrane ; Peripheral membrane protein . Cell junction, tight junction . Recycling endosome . Cell projection . Cytoplasm, cytoskeleton . Cytoplasm, cytosol .
Tissue Specificity	
Function	function:May be a cytoskeletal regulator.,similarity:Contains 1 CH (calponin-homology) domain.,similarity:Contains 1 LIM zinc-binding domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human liver cancer. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).