




ZN503 rabbit pAb

Catalog No	YP-Ab-12168
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZNF503
Immunogen	Synthesized peptide derived from human ZN503 AA range: 469-519
Specificity	This antibody detects endogenous levels of ZN503 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZN503;ZNF503
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:May function as a transcriptional repressor.,similarity:Belongs to the Elbow/Noc family.,similarity:Contains 1 C2H2-type zinc finger.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from PC-12 cells,
primary antibody was diluted at 1:1000, 4° over night