




TM52B rabbit pAb

Catalog No	YP-Ab-12169
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	TMEM52B C12orf59 UNQ5927/PRO19821
Immunogen	Synthesized peptide derived from human TM52B AA range: 128-178
Specificity	This antibody detects endogenous levels of TM52B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	TM52B;TMEM52B C12orf59 UNQ5927/PRO19821
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from VEC cells,
primary antibody was diluted at 1:1000, 4°over night