




ZSWM7 rabbit pAb

Catalog No	YP-Ab-12181
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	ZSWIM7 SWS1
Immunogen	Synthesized peptide derived from human ZSWM7 AA range: 68-118
Specificity	This antibody detects endogenous levels of ZSWM7 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ZSWM7;ZSWIM7 SWS1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	
Function	function:Involved in early stages of the homologous recombination repair (HRR) pathway of double-stranded DNA breaks arising during DNA replication or induced by DNA-damaging agents.,similarity:Belongs to the SWS1 family.,similarity:Contains 1 SWIM-type zinc finger.,subunit:Interacts with RAD51L3.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HpeG2 cells,
primary antibody was diluted at 1:1000, 4°over night