




ARH37 rabbit pAb

Catalog No	YP-Ab-12184
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ARHGEF37
Immunogen	Synthesized peptide derived from human ARH37 AA range: 24-74
Specificity	This antibody detects endogenous levels of ARH37 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ARH37; ARHGEF37
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm,
Tissue Specificity	
Function	function: May act as a guanine nucleotide exchange factor (GEF)., similarity: Contains 1 BAR domain., similarity: Contains 1 DH (DBL-homology) domain., similarity: Contains 2 SH3 domains.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from Hela cells,
primary antibody was diluted at 1:1000, 4° over night