




UD2A3 rabbit pAb

Catalog No	YP-Ab-12201
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	UGT2A3 UNQ2559/PRO6239
Immunogen	Synthesized peptide derived from human UD2A3 AA range: 384-434
Specificity	This antibody detects endogenous levels of UD2A3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UD2A3;UGT2A3 UNQ2559/PRO6239
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Membrane ; Single-pass type I membrane protein .
Tissue Specificity	
Function	catalytic activity:UDP-glucuronate + acceptor = UDP + acceptor beta-D-glucuronoside.,function:UDP-glucuronosyltransferases catalyze phase II biotransformation reactions in which lipophilic substrates are conjugated with glucuronic acid to increase water solubility and enhance excretion. They are of major importance in the conjugation and subsequent elimination of potentially toxic xenobiotics and endogenous compounds.,similarity:Belongs to the UDP-glycosyltransferase family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night