




Z354A rabbit pAb

Catalog No	YP-Ab-12206
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	ZNF354A EZNF HKL1 TCF17
Immunogen	Synthesized peptide derived from human Z354A AA range: 272-322
Specificity	This antibody detects endogenous levels of Z354A at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	EZNF;HKL1;KID 1;KID1;TCF 17;TCF17;Transcription factor 17;zinc finger protein 354A;Zinc finger protein eZNF;ZNF354A;Z354A;ZNF354A EZNF HKL1 TCF17
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in kidney and inner ear.
Function	similarity:Belongs to the krueppel C2H2-type zinc-finger protein family., similarity:Contains 1 KRAB domain.,similarity:Contains 13 C2H2-type zinc fingers.,tissue specificity:Expressed in kidney and inner ear.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night