




GG6L2 rabbit pAb

Catalog No	YP-Ab-12215
Isotype	IgG
Reactivity	Human;Rat;Mouse;
Applications	WB
Gene Name	GOLGA6L2
Immunogen	Synthesized peptide derived from human GG6L2 AA range: 555-605
Specificity	This antibody detects endogenous levels of GG6L2 at Human
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	GG6L2;GOLGA6L2
Observed Band	
Calculated Molecular Weight	
Cell Pathway	atus,
Tissue Specificity	
Function	similarity:Belongs to the GOLGA6 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MDA-MB cells,
primary antibody was diluted at 1:1000, 4° over night