




F190B rabbit pAb

Catalog No	YP-Ab-12229
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM190B KIAA1128 NPD012
Immunogen	Synthesized peptide derived from human F190B AA range: 609-659
Specificity	This antibody detects endogenous levels of F190B at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1128;Gcap14;Protein FAM190B;Protein GCAP14 homolog;F190B; FAM190B KIAA1128 NPD012
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytoskeleton . Associates with microtubules in interphase. Has diffuse expression throughout the cell during mitosis. .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night