




P20D1 rabbit pAb

Catalog No	YP-Ab-12230
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	PM20D1
Immunogen	Synthesized peptide derived from human P20D1 AA range: 3-53
Specificity	This antibody detects endogenous levels of P20D1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Peptidase M2 domain-containing protein 1;N-fatty-acyl-amino acid synthase/hydrolase PM2D1;EC:3.5.1.114;P20D1;PM20D1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	cofactor: Binds 2 zinc ions per subunit.,similarity: Belongs to the peptidase M20A family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from 293T cells,
primary antibody was diluted at 1:1000, 4° over night