




FA13C rabbit pAb

Catalog No	YP-Ab-12232
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	FAM13C FAM13C1
Immunogen	Synthesized peptide derived from human FA13C AA range: 422-472
Specificity	This antibody detects endogenous levels of FA13C at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FA13C;FAM13C FAM13C1
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	similarity:Belongs to the FAM13 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from AD293 cells,
primary antibody was diluted at 1:1000, 4° over night