




DMRT3 rabbit pAb

Catalog No	YP-Ab-12246
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	DMRT3 DMRTA3
Immunogen	Synthesized peptide derived from human DMRT3 AA range: 219-269
Specificity	This antibody detects endogenous levels of DMRT3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Doublesex- and mab-3-related transcription factor 3;DMRTA3;DMRT3;DMRT3 DMRTA3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Nucleus .
Tissue Specificity	Expressed in testis.
Function	developmental stage:Expressed in 4 to 5 weeks embryos.,function:May regulate transcription during sexual development.,similarity:Belongs to the DMRT family., similarity:Contains 1 DM DNA-binding domain.,tissue specificity:Specifically expressed in testis.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from HuvEc cells,
primary antibody was diluted at 1:1000, 4°over night