



JOS2 rabbit pAb


Catalog No	YP-Ab-12248
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB;IHC
Gene Name	JOSD2 SBBI54
Immunogen	Synthesized peptide derived from human JOS2 AA range: 62-112
Specificity	This antibody detects endogenous levels of JOS2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	FLJ29018;JOS2_HUMAN;JOSD 2;JOSD2;Josephin 2;Josephin domain containing 2;Josephin domain containing protein 2;Josephin domain-containing protein 2;Josephin-2;Josephin2;SBBI54;JOS2;JOSD2 SBBI54
Observed Band	21kD
Calculated Molecular Weight	21kD
Cell Pathway	Cytoplasm, cytosol .
Tissue Specificity	
Function	similarity:Contains 1 Josephin domain.,
Background	This gene encodes a protein containing a Josephin domain. Josephin domain-containing proteins are deubiquitinating enzymes which catalyze the hydrolysis of the bond between the C-terminal glycine of the ubiquitin peptide and protein substrates. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Jul 2012],
matters needing attention	Avoid repeated freezing and thawing!




Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from DU145 cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human Gastric adenocarcinoma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).