




CA177 rabbit pAb

Catalog No	YP-Ab-12259
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB
Gene Name	C1orf177
Immunogen	Synthesized peptide derived from human CA177 AA range: 300-350
Specificity	This antibody detects endogenous levels of CA177 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	CA177;C1orf177
Observed Band	
Calculated Molecular Weight	
Cell Pathway	mitochondrion,mitochondrial ribosome,
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from MCF-7 cells,
primary antibody was diluted at 1:1000, 4° over night