



MRI rabbit pAb


Catalog No	YP-Ab-12270
Isotype	IgG
Reactivity	Human; Mouse;Rat
Applications	WB;ELISA;IHC
Gene Name	MRI C7orf49
Immunogen	Synthesized peptide derived from human MRI AA range: 40-90
Specificity	This antibody detects endogenous levels of MRI at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000;IHC-p 1:50-300; ELISA 2000-20000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	Cell cycle regulator of NHEJ;Cell cycle regulator of non-homologous end joining; CYREN;MRI;MRI C7orf49
Calculated Molecular Weight	17kD
Cell Pathway	[Isoform 1]: Cytoplasm . Nucleus . Chromosome . Nuclear localization may depend upon interaction with XRCC5/Ku80 and XRCC6/Ku70 heterodimer (PubMed:24610814). Localizes to DNA damage sites (PubMed:27063109). . ; [Isoform 3]: Cytoplasm . Some nuclear localization may be due to passive diffusion. . ; [Isoform 4]: Cytoplasm . Nucleus . Nuclear localization may depend upon interaction with XRCC5/Ku80 and XRCC6/Ku70 heterodimer and increases upon etoposide treatment. .
Tissue Specificity	
Function	function:May act as a regulator of proteasome.,
Background	
matters needing attention	Avoid repeated freezing and thawing!



Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images



Western blot analysis of lysates from 293T cells, primary antibody was diluted at 1:1000, 4° over night



Immunohistochemical analysis of paraffin-embedded human oophoroma. 1, Antibody was diluted at 1:200(4° overnight). 2, Tris-EDTA, pH9.0 was used for antigen retrieval. 3, Secondary antibody was diluted at 1:200(room temperature, 45min).