




CC157 rabbit pAb

Catalog No	YP-Ab-12281
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	CCDC157 KIAA1656
Immunogen	Synthesized peptide derived from human CC157 AA range: 16-66
Specificity	This antibody detects endogenous levels of CC157 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	KIAA1656;Coiled-coil domain-containing protein 157;CC157;CCDC157 KIAA1656
Observed Band	
Calculated Molecular Weight	
Cell Pathway	
Tissue Specificity	
Function	sequence caution:Intron retention. This sequence is incomplete at the 5'- and 3'-ends and extensively differs from that shown due to multiple introns.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from U2OS cells,
primary antibody was diluted at 1:1000, 4° over night