




ORML3 rabbit pAb

Catalog No	YP-Ab-12283
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	ORMDL3
Immunogen	Synthesized peptide derived from human ORML3 AA range: 98-148
Specificity	This antibody detects endogenous levels of ORML3 at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	ORM1 (S. cerevisiae)-like 3; ORM1-like 3 (S. cerevisiae); ORM1-like 3; ORM1-like protein 3; ORMDL sphingolipid biosynthesis regulator 3; Ormdl3; ORML3_HUMAN; OTTHUMP00000164334; ORML3
Observed Band	17kD
Calculated Molecular Weight	17kD
Cell Pathway	Endoplasmic reticulum membrane ; Multi-pass membrane protein .
Tissue Specificity	Widely expressed. Expressed in adult and fetal heart, brain, lung, liver, skeletal muscle and kidney. Expressed in adult pancreas and placenta and in fetal spleen and thymus.
Function	similarity: Belongs to the ORM family., tissue specificity: Widely expressed. Expressed in adult and fetal heart, brain, lung, liver, skeletal muscle and kidney. Expressed in adult pancreas and placenta and in fetal spleen and thymus.,
Background	
matters needing attention	Avoid repeated freezing and thawing!

Usage suggestions

This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.

Products Images

Western blot analysis of lysates from CACO2 cells, primary antibody was diluted at 1:1000, 4° over night