




BPIA3 rabbit pAb

Catalog No	YP-Ab-12291
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	BPIFA3 C20orf71 SPLUNC3
Immunogen	Synthesized peptide derived from human BPIA3 AA range: 73-123
Specificity	This antibody detects endogenous levels of BPIA3 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	BPIA3;BPIFA3 C20orf71 SPLUNC3
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	similarity:Belongs to the BPI/LBP/Plunc superfamily. Plunc family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A549 cells,
primary antibody was diluted at 1:1000, 4° over night