




NXNL2 rabbit pAb

Catalog No	YP-Ab-12294
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	NXNL2 C9orf121
Immunogen	Synthesized peptide derived from human NXNL2 AA range: 36-86
Specificity	This antibody detects endogenous levels of NXNL2 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	C9orf121;Chromosome 9 open reading frame 121;Nucleoredoxin like 2; Nucleoredoxin-like protein 2;Nxnl2;NXNL2_HUMAN;RDCVF2;NXNL2 C9orf121
Observed Band	
Calculated Molecular Weight	
Cell Pathway	cytoplasm,
Tissue Specificity	
Function	similarity:Belongs to the nucleoredoxin family.,similarity:Contains 1 thioredoxin domain.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from HeLa cells,
primary antibody was diluted at 1:1000, 4° over night