




UTS2B rabbit pAb

Catalog No	YP-Ab-12295
Isotype	IgG
Reactivity	Human; Mouse; Rat
Applications	WB
Gene Name	UTS2D URP UTS2B
Immunogen	Synthesized peptide derived from human UTS2B AA range: 53-103
Specificity	This antibody detects endogenous levels of UTS2B at Human/Mouse/Rat
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit, IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	UTS2B; UTS2D URP UTS2B
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Secreted .
Tissue Specificity	
Function	function: Potent vasoconstrictor., similarity: Belongs to the urotensin-2 family.,
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from A431 cells,
primary antibody was diluted at 1:1000, 4° over night