




AP5S1 rabbit pAb

Catalog No	YP-Ab-12303
Isotype	IgG
Reactivity	Human; Mouse
Applications	WB
Gene Name	AP5S1 C20orf29
Immunogen	Synthesized peptide derived from human AP5S1 AA range: 146-196
Specificity	This antibody detects endogenous levels of AP5S1 at Human/Mouse
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Source	Polyclonal, Rabbit,IgG
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution	WB 1:500-2000
Concentration	1 mg/ml
Purity	≥90%
Storage Stability	-20°C/1 year
Synonyms	AP5S1;AP5S1 C20orf29
Observed Band	
Calculated Molecular Weight	
Cell Pathway	Cytoplasm, cytosol . Late endosome membrane ; Peripheral membrane protein ; Cytoplasmic side . Lysosome membrane ; Peripheral membrane protein ; Cytoplasmic side .
Tissue Specificity	
Function	
Background	
matters needing attention	Avoid repeated freezing and thawing!
Usage suggestions	This product can be used in immunological reaction related experiments. For more information, please consult technical personnel.



Products Images



Western blot analysis of lysates from DU145 cells,
primary antibody was diluted at 1:1000, 4° over night